

# -

## Product Information Sheet

## **FITC Conjugated Protein A**

Catalog No.	BA1120-1
outurog no.	D/(1120 1

1 mg

#### Storage

Size

Store at 4 °C. Avoid light and multiple freeze-thaw cycles.

#### Expiration

One year from the day of shipment.

#### Applications

Flowcytometry (FCM); Immunocytochemistry (IC); Indirect Immunofluorescence (IF); Immunohistochemistry (IHC).

#### **Product Description**

Protein A is the component of the cytoderm of staphylococcus aureus. It was regarded as a kind of broad spectrum secondary antibody, because it can bind Fc fragments of a majority of immunoglobulins. It can react with guinea pig, rabbit globulin and human IgG1, 2, 4 strongly and react with mouse, bovine globulin and human IgA,E,M weakly.

#### Contents

1 mg of FITC conjugated protein A, F/P=3/1; 0.01M PBS; 1% BSA; 0.01% Thimerosal.

#### Labeling Method

Protein A is conjugated to FITC by means of a method described by Hijmans,W.,et al.

(Reference: Hijmans, W., et al. Clin. Exp. Immunol., 4, 457(1969))

#### Preparation of Diluent Buffer

Add reagent grade BSA into 0.01 M PBS (PH7.2-7.6) or TBS buffer and make BSA at a concentration of 1%. Use the above diluent buffer to dilute. See "Recommended Dilutions" below for details.

Preparation of 0.01M **TBS:** Add 1.2g Tris, 8.5g NaCl;  $450\mu$ I of purified acetic acid or 700 $\mu$ I of concentrated hydrochloric acid to 1000ml H<sub>2</sub>O and adjust pH to 7.2-7.6. Finally, adjust the total volume to 1L.

Preparation of 0.01M **PBS:** Add 8.5g sodium chloride,  $1.4g Na_2HPO_4$  and 0.2g NaH<sub>2</sub>PO<sub>4</sub> to 1000ml distilled water and adjust pH to 7.2-7.6. Finally, adjust the total volume to 1L.

Recommended Dilutions

To reorder contact us at:	FC	
Antagene, Inc.	IC	
I Free: 1(866)964-2589	Inc	
Tel: (650) 964-2589	IH	
Fax: (650) 964-2519	Οµ	
nail: Info@antageneinc.com		

FCM	1:100
IC	1:64-128
Indirect IF	1:64-128
IHC (F)	1:64-128

Optimal working dilutions must be determined by end user.

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC AND CLINICAL USE.

Tol

em

# SOLID CONVERTER