



## **Product Information Sheet**

## 10 nm Colloidal Gold Conjugated Goat Anti-mouse IgG

Catalog No. GA1004

Purification This antibody was purified from antiserum by immunoaffinity chromatography which

Size 0.1 mg

removes essentially all goat serum proteins, except the specific antibody for mouse IgG.

Ig Type IgG

for at least one year.

Raised in Goat Clonality Polyclonal

Immunogen

Mouse IgG (whole molecular).

**Contents** 

0.1 mg of colloidal gold conjugated specific antibody, the diameter of colloidal gold is Storage

10±2 nm; 0.01M PBS; 1% BSA; 0.01% Thimerosal.

Avoid multiple freeze-thaw cycles.

Specificity

This colloidal gold conjugated antibody is specific for mouse IgG and shows no **Expiration** 

cross-reactivity with mouse IgA, IgM and other animal species IgG proteins.

**Applications** 

One year from the day of shipment.

Immunohistochemistry (IHC).

Preparation of Colloidal Gold

10 nm Colloidal Gold - Add 0.05% trisodium citrate into boiled 0.01% chloroauric acid Electron Microscopy (EM);

which was reduced to gold granule solution through redox reaction.

Conjugate colloidal gold to antibody after mensurating the optimal data. And through overspeed centrifugation and gradient delamination to density, the colloidal gold

conjugated reagents were produced.

Preparation of Diluent Buffer

Add reagent grade BSA into 0.01 M PBS (PH7.2-7.6) or TBS buffer and make BSA at a concentration of 1%. Use the above diluent buffer to dilute. See "Recommended

Dilutions" below for details.

Preparation of IGSS (ImmunoGold Silver Staining)

Add citric acid 2.35g, trisodium citrate 2.55g, glutin 1g and hydroquinone 1.7g into 100 ml of distilled water and mix evenly. Before usage, add silver nitrate 50mg into it. Then,

add this mixed solution into tissue at last and colorate for 10-20 minutes in dark. Wash

with water to terminate reaction.

To reorder contact us at: Antagene, Inc.

Toll Free: 1(866)964-2589

Tel: (650) 964-2589

Fax: (650) 964-2519

email: Info@antageneinc.com

Recommended Dilutions

ΕM 1:20-50

IHC 1:50-100

Optimal working dilutions must be determined by end user.