



Product Information Sheet

FITC Conjugated Goat Anti-rabbit IgG

Catalog No. BA1105-1

Size 1 mg

Ig Type IgG

Immunogen

Rabbit IgG (whole molecular).

Storage

Store at 4 °C. Avoid light and multiple freeze-thaw cycles.

Expiration

One year from the day of shipment.

Applications

Flowcytometry (FCM);

Immunocytochemistry (IC);

Indirect Immunofluorescence (IF); Immunohistochemistry (IHC). **Purification**

This antibody was purified from antiserum by immunoaffinity chromatography which removes essentially all goat serum proteins, except the specific antibody

for rabbit IgG.

Raised in Goat

Clonality Polyclonal

Contents

1 mg of FITC conjugated specific antibody (purity is above 99%), FITC : Ab = 4-6:1;0.01M PBS; 1% BSA; 0.01% Thimerosal. The emission and filtration

wavelength of FITC are 495 nm and 525 nm respectively.

Specificity

This FITC conjugated antibody is specific for rabbit IgG and shows no cross-reactivity with rabbit IgA, IgM and other animal species IgG proteins.

Labeling Method

Goat anti-rabbit IgG is conjugated to FITC by means of a method described by Hijmans,W.,et al.

(Reference: Hijmans, W., et al. Clin. Exp. Immunol., 4, 457(1969))

Preparation of Diluent Buffer

Add reagent grade BSA into 0.01 M PBS (PH7.2-7.6) or TBS buffer and make BSA at a concentration of 1%. Use the above diluent buffer to dilute. See "Recommended Dilutions" below for details.

Preparation of 0.01M **TBS:** Add 1.2g Tris, 8.5g NaCl; 450 μ l of purified acetic acid or 700 μ l of concentrated hydrochloric acid to 1000ml H₂O and adjust pH to 7.2-7.6. Finally, adjust the total volume to 1L.

Preparation of 0.01M **PBS:** Add 8.5g sodium chloride, 1.4g Na₂HPO₄ and 0.2g NaH₂PO₄ to 1000ml distilled water and adjust pH to 7.2-7.6. Finally, adjust the total volume to 1L.

c.com Recommended Dilutions

FCM 1:100

To reorder contact us at:

Antagene, Inc.

Toll Free: 1(866)964-2589

Tel: (650) 964-2589

Fax: (650) 964-2519

email: Info@antageneinc.com

Product Information Sheet

IC	1:32-64	IHC	1:32-64
Indirect IF	1:32-64	Optimal working dilutions must be determined by end user.	